



## LOIRE Toilet Support Arms

Oval shaped profile is great to push on and pull against  
Easy to raise and lower  
Very stable in use for extra confidence  
Easy to clean

- The wave design or the Raised Attachments can virtually eliminate the need for height adjustment.
- Each has high level support at the front, which when set at 80cm high, meets the needs of nearly all standing transfers. This high support is clear of the front of the toilet pan.
- The lower levels alongside the toilet provide a second level for easy lowering or standing, while providing a good height for side transfers.
- When sitting, the lower height makes dressing and undressing easier.

# LOIRE Toilet Support Arms

LOIRE toilet support arms in two versions - Straight or Wave profiles.  
Both Straight and Wave are delivered in length 90 cm

<p>LOIRE-straight 80 Item 40-40110</p>  <p>Max load: 80 kg Length: 90 cm</p>	<p>LOIRE-wave 80 Item 40-40115</p>  <p>Max load: 80 kg Length: 90 cm</p>	<p>Get up/sit down attachment Item 40-40925</p> 
<p>LOIRE-straight 135, with leg Item 40-40120, telescopic leg</p>  <p>Max load: 135 kg Length: 90 cm</p>	<p>LOIRE-wave 135, with leg Item 40-40125, telescopic leg</p>  <p>Max load: 135 kg Length: 90 cm</p>	<p>Towel holder Item 40-40915</p> 
<p>LOIRE-straight 140 Item 40-40130</p>  <p>Max load: 140 kg Length: 90 cm</p>	<p>LOIRE-wave 140 Item 40-40135</p>  <p>Max load: 140 kg Length: 90 cm</p>	<p>Toilet roll holder Item 40-40910</p> 
<p>Wall Plate for height adjustment: LOIRE: Item 40-40905</p> 		

## ROPOX Tips

- Fix front of rail at 80cm from floor for most standing transfers.
- If using straight rail with Get Up/Sit Down attachment, place rail at 68cm from floor.
- For sideways pulling on rails, the Loire offers great stability.
- If height adjustment is still considered necessary, the Wall Plate accessory can be used with all the models to give 12cm of adjustment.